



Official Website



JUMPER Health

Shenzhen Jumper Medical Equipment Co., Ltd.

Address: D Building, No. 71, Xintian Road, Fuyong Street, Baoan, Shenzhen, Guangdong, China

Tel: +86 755-26696279 E-mail: info@jumper-medical.com Web: www.jumper-medical.com

**PRODUCT
BROCHURE**

**SINCE
1998**



JUMPER

Healthcare within mind

www.jumper-medical.com



BLOOD PRESSURE MONITOR



-  Bluetooth
-  Guest Mode
-  2 Users
-  Voice Switch
-  Type-C
-  199 Records
-  Average Value
-  30s Auto Power Off
-  Irregular Heartbeat Indicator

FIND THE RIGHT BLOOD PRESSURE MONITOR FOR YOU.

Our blood pressure monitor features advanced technology to automatically inflate and deflate at the appropriate level, delivering a personalized accurate, and comfortable measurement.

It can also calculate average blood pressure value based on the last three readings. It is convenient for you to use at home or take it out because of the compact size, so you'll never miss a measurement.

FEATURES

- Universal cuff
- Extra-large screen
- One touch operation
- Display date and time
- Lightweight & portable
- Average of last 3 readings
- Irregular heartbeat indicator
- Auto-off power saving feature
- 2-User mode & memory function





HAA11

BLUETOOTH	WITHOUT BLUETOOTH	TYPE	LCD SIZE
HAA11	HAA10	ARM TYPE	91 x 77 mm / 4.8 inch

USER MEMORY	AUTOMATIC SHUTDOWN	CUFF LENGTH(cm)
2 x 199 RECORDS	30s	22-36 / 22-42

MEASUREMENT RANGE	POWER SUPPLY
PRESSURE VALUE 0-295 ± 3 mmHg PULSE RATE 40-199 ± 5% BPM	AA x 4 USB 5V



HAA21

BLUETOOTH	WITHOUT BLUETOOTH	TYPE	LCD SIZE
HAA21	HAA20	ARM TYPE	71 x 101 mm / 4.9 inch

USER MEMORY	AUTOMATIC SHUTDOWN	CUFF LENGTH(cm)
2 x 199 RECORDS	30s	22-36 / 22-42

MEASUREMENT RANGE	POWER SUPPLY
PRESSURE VALUE 0-295 ± 3 mmHg PULSE RATE 40-199 ± 5% BPM	AAA x 4 USB 5V




JPD-HA101 (LCD)

TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	83 x 97 mm / 5.1 inch	22-36 / 22-42

USER MEMORY	AUTOMATIC SHUTDOWN
2 x 99 RECORDS	30s

MEASUREMENT RANGE	POWER SUPPLY
PRESSURE VALUE 0-295 ± 3 mmHg PULSE RATE 40-199 ± 5% BPM	AA x 4 USB 5V



JPD-HA101 (LED)

TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	83 x 97 mm / 5.1 inch	22-36 / 22-42

USER MEMORY	AUTOMATIC SHUTDOWN
2 x 99 RECORDS	30s

MEASUREMENT RANGE	POWER SUPPLY
PRESSURE VALUE 0-295 ± 3 mmHg PULSE RATE 40-199 ± 5% BPM	AA x 4 USB 5V




JPD-HA121

TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	83 x 97 mm / 5.1 inch	22-36 / 22-42

USER MEMORY	AUTOMATIC SHUTDOWN
2 x 99 RECORDS	30s

MEASUREMENT RANGE	POWER SUPPLY
PRESSURE VALUE 0-295 ± 3 mmHg PULSE RATE 40-199 ± 5% BPM	AA x 4 USB 5V



A

JPD-HA100


TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	60 x 77 mm / 3.9 inch	22-36 / 22-42
USER MEMORY	AUTOMATIC SHUTDOWN	
2 x 99 RECORDS	30s	
MEASUREMENT RANGE		POWER SUPPLY
PRESSURE VALUE	0-295 ± 3 mmHg	AA x 4
PULSE RATE	40-199 ± 5% BPM	USB 5V



JPD-HA200




TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	59 x 47 mm / 3 inch	22-36 / 22-42
USER MEMORY	AUTOMATIC SHUTDOWN	
2 x 99 RECORDS	30s	
MEASUREMENT RANGE		POWER SUPPLY
PRESSURE VALUE	0-295 ± 3 mmHg	AAA x 4
PULSE RATE	40-199 ± 5% BPM	USB 5V

**B**

JPD-HA100

TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	60 x 77 mm / 3.9 inch	22-36 / 22-42
USER MEMORY	AUTOMATIC SHUTDOWN	
2 x 99 RECORDS	30s	
MEASUREMENT RANGE		POWER SUPPLY
PRESSURE VALUE	0-295 ± 3 mmHg	AA x 4
PULSE RATE	40-199 ± 5% BPM	USB 5V



JPD-HA210



TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	59 x 47 mm / 3 inch	22-36 / 22-42
USER MEMORY	AUTOMATIC SHUTDOWN	
2 x 99 RECORDS	30s	
MEASUREMENT RANGE		POWER SUPPLY
PRESSURE VALUE	0-295 ± 3 mmHg	AAA x 4
PULSE RATE	40-199 ± 5% BPM	USB 5V



JPD-HA120




TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	60 x 77 mm / 3.9 inch	22-36 / 22-42
USER MEMORY	AUTOMATIC SHUTDOWN	
2 x 99 RECORDS	30s	
MEASUREMENT RANGE		POWER SUPPLY
PRESSURE VALUE	0-295 ± 3 mmHg	AA x 4
PULSE RATE	40-199 ± 5% BPM	USB 5V



JPD-HA300



TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	59 x 47 mm / 3 inch	22-36 / 22-42
USER MEMORY	AUTOMATIC SHUTDOWN	
2 x 99 RECORDS	30s	
MEASUREMENT RANGE		POWER SUPPLY
PRESSURE VALUE	0-295 ± 3 mmHg	AAA x 4
PULSE RATE	40-199 ± 5% BPM	USB 5V



NEW ARRIVAL




HWA10

TYPE: WRIST
LCD SIZE: 45 x 41 mm / 2.4 inch
CUFF LENGTH(cm): 13.5-21.5

USER MEMORY: 2 x 199 RECORDS
POWER SUPPLY: AAA x 2
AUTO SHUTDOWN: 30s

MEASUREMENT RANGE
PRESSURE VALUE: 0-295 ± 3 mmHg
PULSE RATE: 40-199 ± 5% BPM



NEW ARRIVAL




HWA20

TYPE: WRIST
LCD SIZE: 41x 41 mm / 2.3 inch
CUFF LENGTH(cm): 13.5-21.5

USER MEMORY: 2 x 199 RECORDS
POWER SUPPLY: AAA x 2
AUTO SHUTDOWN: 30s

MEASUREMENT RANGE
PRESSURE VALUE: 0-295 ± 3 mmHg
PULSE RATE: 40-199 ± 5% BPM



NEW ARRIVAL




HWA11

TYPE: WRIST
LCD SIZE: 45 x 41 mm / 2.4 inch
CUFF LENGTH(cm): 13.5-21.5

USER MEMORY: 2 x 199 RECORDS
POWER SUPPLY: AAA x 2
AUTO SHUTDOWN: 30s

MEASUREMENT RANGE
PRESSURE VALUE: 0-295 ± 3 mmHg
PULSE RATE: 40-199 ± 5% BPM



NEW ARRIVAL




HWA21

TYPE: WRIST
LCD SIZE: 41x 41 mm / 2.3 inch
CUFF LENGTH(cm): 13.5-21.5

USER MEMORY: 2 x 199 RECORDS
POWER SUPPLY: AAA x 2
AUTO SHUTDOWN: 30s

MEASUREMENT RANGE
PRESSURE VALUE: 0-295 ± 3 mmHg
PULSE RATE: 40-199 ± 5% BPM



JPD-900W



TYPE: WRIST TYPE
LCD SIZE: 37 x 45 mm / 2.3 inch
CUFF LENGTH(cm): 13.5-19.5

USER MEMORY: 90 RECORDS
AUTOMATIC SHUTDOWN: 1 min

MEASUREMENT RANGE
PRESSURE VALUE: 30-280 ± 3 mmHg
PULSE RATE: 40-199 ± 5% BPM
POWER SUPPLY: LR03 AAA x 2










Blood Pressure Monitor

Parameters



Product Model	HAA11	HAA21	HAA10	HAA20	JPD-HA300	JPD-HA210		JPD-HA200	JPD-HA120	JPD-HA121	JPD-HA100	JPD-HA100	JPD-HA101(LCD)	JPD-HA101(LED)
Type	Arm							Arm						
Pressure Value	0-295 ± 3 mmHg							0-295 ± 3 mmHg						
Pulse Rate	40-199 ± 5% BPM							40-199 ± 5% BPM						
Display Size	mm	91 x 77	71 x 101	91 x 77	71 x 101	59 x 47		59 x 47	60 x 77	83 x 97	60 x 77		83 x 97	
	inch	4.8	4.9	4.8	4.9	3		3	3.9	5.1	3.9		5.1	
Average	✓							✓						
Bluetooth	✘		-				-	✘		-				
Backlight	☀			-			-	☀		-		☀		-
Auto Power Off	30s							30s						
Cuff Length (cm)	22-36 / 22-42 (Optional)							22-36 / 22-42 (Optional)		22-36 / 22-42 (Optional)				
Power Supply	AA x 4 USB 5V	AAA x 4 USB 5V	AA x 4 USB 5V	AAA x 4 USB 5V				AAA x 4 USB 5V	AA x 4 USB 5V					
User Memory	2 x 199 records				2 x 99 records			2 x 99 records	2 x 99 records					
Type-C	☐				-			-						☐

Blood Pressure Monitor							
Parameters							
 Product Model	JPD-900A	HWA11	HWA21	HWA10	HWA20	JPD-900W	
 Type	Arm	Wrist					
 Pressure Value	30-280 ± 3 mmHg	0-295 ± 3 mmHg		0-295 ± 3 mmHg		30-280 ± 3 mmHg	
 Pulse Rate	40-199 ± 5% BPM						
 Display Size	mm	50 × 60	45 × 41	41 × 41	45 × 41	41 × 41	37 × 45
	inch	3.1	2.4	2.3	2.4	2.3	2.29
 Average	—	✓				—	
 Bluetooth	—	✕		—			
 Backlight	—	☀				—	
 Auto Power Off	1 min	30s				1 min	
 Cuff Length (cm)	22-32	13.5-21.5				13.5-19.5	
 Power Supply	AA x 4 LR6	AAA x 2				AAA x 2 LR03	
 User Memory	90 records	2 x 199 records				90 records	
 Type-C	—						

